Description

Light-curing, single-component bonding resin for optimizing the enamel-etch technique in combination with all light-curing restorative materials.

Composition

Heliobond contains:

Bis-GMA	60 %wt.
---------	---------

Triethylene glycol dimethacrylate 40 %wt.

Indication

- Bonding agent for adhesive restorations.
- Transparent sealant for fissures and restorations.
- Bonding layer for the repair of resin-, crown-, or bridgework.

Contraindication

Heliobond should not be used if patients are known to be allergic to any of the ingredients.

Interactions

To achieve an optimum bond to the enamel, the enamel-etch technique must be applied. Contamination of the etched and dried enamel surface with saliva or blood will hamper this effect. Phenolic substances (e.g. eugenol) inhibit polymerization. The use of base materials containing these substances must be avoided.

Application

Preparation of enamel

- Depending on the case, bevel enamel surfaces to be treated using a diamond or clean them carefully with an abrasive (e.g. pumice powder and water). Do not use any oily or greasy paste. Possible residue from base materials on the enamel should be removed in this way. Then rinse and dry.
- Apply Total Etch on the enamel surfaces to be treated and allow it to react for 15–30 s. Afterwards, carefully rinse with water spray and dry the etched surfaces with a stream of air (use only water- and oil-free air).
- A dry operating field must be ensured. Ideally, a rubber dam should be placed.
- The etched and dried enamel surface must not be contaminated prior to the application of Heliobond (if the surface comes into contact with saliva or blood, the etching and drying procedures must be repeated).

Dosage

- Dosage depends on the indication. When used as a bonding agent, a very thin layer of Heliobond should be applied. If necessary, use a stream of air to achieve an optimal thin layer.
- When used as a transparent fissure sealant, Heliobond must be applied using a thicker layer.

Heliobond

Heliobond

EN Instructions for Use	– Light-curing bonding resin
DE Gebrauchsinformation	– Lichthärtender Haftvermittler
FR Mode d'emploi	– Agent de liaison photopolymérisable
IT Istruzioni d'uso	– Promotore di adesione fotoindurente
ES Instrucciones de uso	– Adhesivo fotopolimerizable
PT Instruções de Uso	– Adesivo resinoso fotopolimerizável
SV Bruksanvisning	– Ljushärdande bonding resin
DA Brugsanvisning	– Lys hærdende resin
FI Käyttöohjeet	– Valokovettinen killesidosaine
NO Bruksanvisning	– Lys herdende bonding
NL Gebruiksaanwijzing	– Lichtuithardend hechtmiddel
EL Οδηγίες Χρήσεως	– Φωτοπολυμεριζούμενη υγρή ρητίνη
TR Kullanma Talimatı	– Işıklı sertleşen bonding ajanı
SL Navodila za uporabo	– Pomožno adhezivno sredstvo, strjuje se na svetlobi
HR Upute za uporabu	– Svjetlom polimerizirajuće vezivno sredstvo
CS Návod k použití	– Světlem tuhnoucí adhezivum
SK Návod na používanie	– Svetlom vytvrdzovaný bondovací prípravok
HU Használati utasítás	– Fényre keményedő kötőanyagta
PL Instrukcja stosowania	– Światłoutwardzalny, jednoskładnikowy materiał łączący
SR Упутство за употребу	– Светлосно полимеризујућа везивна смола
MK Упатство за употреба	– Светлосно-полимеризирачка смола за бондирање
BG Инструкции за употреба	– Фотополимеризиращо свързващо средство
SQ Udhëzime përdorimi	– Preparaat adeziv fotopolimerizues
RO Instrucțiuni de utilizare	– Agent de legătură fotopolimerizabil
RU Инструкция по применению	– Светоотверждаемый композиционный бондинг
UK Інструкції до застосування	– Адгезивна смола
ET Kasutamisi juhend	– Valguskõvastuv vaik
LV Lietošanas instrukcija	– Gaismā cietējoša sveķu saistviela
LT Naudojimo instrukcija	– Šviesoje kietėjanti rišamoji derva

Method of application

1. Bonding agent for adhesive restorations

- Etch desired enamel surfaces (see preparation of the enamel). If necessary, cover exposed dentin surfaces with Syntac® Adhesive (see Syntac instructions for use) or appropriate base material.
- Apply a thin layer of Heliobond onto the etched enamel surface using a brush or spherical instrument.
- An optimal, thin layer can be achieved using a stream of air.
- If used in combination with light-curing luting composites, Heliobond does not require separate polymerization. For direct restorations, separate polymerization of Heliobond is required: Light-cure for 10 s using a light intensity of 500 to 1100 mW/cm².
- Application of the composite.
- Polymerization of the composite.
- Finishing of the restoration.

2. Transparent sealant for fissures and pits

- Etch desired enamel surfaces (see preparation of the enamel).
- Apply Heliobond into the fissures using a suitable instrument or brush. Avoid entrapment of air and wait for approx. 15 s to allow penetration.
- Polymerize with light for 20 s using a light intensity of 500 to 1100 mW/cm².
- After curing, remove inhibited layer. Check occlusion and grind off any excess material.

3. Preparation of resin repair work

- Roughen the smooth resin surfaces.
- If necessary, apply Monobond® Plus (see Monobond Plus instructions for use).
- Apply a thin layer of Heliobond (use a stream of air, if necessary).
- Polymerize with light for 10 s using a light intensity of 500 to 1100 mW/cm².

Polymerization

- A polymerization time of 10 s is sufficient for thin layers using a light intensity of 500 to 1100 mW/cm².
- Thicker layers (e.g. sealings) require 20 s curing time using a light intensity of 500 to 1100 mW/cm².
- The light probe should be held less than 5 mm above the Heliobond surface and should not touch unpolymerized material.

Note

- Do not expose Heliobond to intensive light during application since this will significantly shorten the working time.
- If Heliobond is used as a cover layer (e.g. sealant), a thin layer of unpolymerized, sticky material will remain on the surface after polymerization due to the oxygen inhibition. It is recommended to remove this layer with cotton rolls, pellets or by polishing.

Storage

- Close bottle immediately after use. Exposure to light may lead to premature polymerization.
- Storage: 2–28 °C / 36–82 °F.
- Do not use Heliobond after the date of expiration.
- Date of expiration: See note on the bottle/packaging.

Warning

Avoid contact of Heliobond with the skin/mucous membrane or eyes. Unpolymerized Heliobond may cause slight, reversible irritation, or lead to a general sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Keep material out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.