OptraGate® Dry Control

[en] Instructions for Use

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1 Intended use

Intended purpose

Suction of liquids in the oral cavity for moisture control, as well as tongue retraction during dental treatment.

Patient target group

- Patients with permanent teeth

Intended users

- Dentists
- Dental assistants
- Dental hygienists / Prophylaxis assistants

Special training

None

Usage


For dental use only.

Description

OptraGate® Dry Control is a dental saliva and liquid ejector with integrated tongue retractor for efficient moisture control in the oral cavity.

The product consists of the following components:

- (1) Tongue retractor (2) Folding hinge, (3) Suction element, (4) Flexible hose (6.3 mm in diameter)

Indications

None

Areas of application:

Moisture control and tongue retraction, e.g. during the following dental interventions:

- Dental examinations
- Partial intraoral scanning
- Tooth preparation
- Direct restorative treatments
- Placing of indirect restorations

Contraindications

- The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Restrictions of use

The product must not be used in the following cases:

- Oral surgical and implantological procedures in the oral cavity
- Sterile procedures
- Reprocessing the product

Do not re-use

Side effects

None known to date

Interactions

None known to date

Clinical benefit

- Moisture control in the treatment area
- Retraction of the patient's tongue

Composition

Polyvinyl chloride (PVC), thermoplastic elastomer (TPE), steel 1.4310, polypropylene (PP)

2 Application

Preparing the product

1. Open the pouch at the marked location.
2. Remove the product from the pouch.
3. Connect the product to the extraction unit or the suitable interface on the treatment unit.

Positioning the product in the oral cavity

1. In very small oral cavities: If necessary, shorten the suction element (3) by approximately 5 mm at the designated position (see image)
2. If required: Position the lip / cheek retractor in the patient's mouth.
3. Position the suction element (3) in the lower jaw next to the tongue.
4. Secure the product to the lower jaw by means of the flexible hose (4).
5. Bend the folding hinge (2) of the tongue retractor (1) according to the size of the palate and the height of the oral cavity.

Removing the product after treatment

1. Remove the product from the oral cavity after treatment.
2. Dispose of the contaminated product.

3 Safety information

- The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- Cleaning and disinfecting agents may negatively affect the material and lead to undesired reactions in patients.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

Supporting documents

Document	To be found at:
Current version of Instructions for Use	www.ivoclar.com/elFU
Structure of Instructions for use and Warnings	www.ivoclar.com/elFU
Explanation of symbols	www.ivoclar.com/elFU

Disposal information

Dispose of the contaminated product according to the respective national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:
None known to date

4 Shelf life and storage

- Protect the product from direct sunlight.
- Date of expiration: See note on packaging

5 Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from improper use of the product. Furthermore, the user is responsible for determining the product's suitability for the user's own purposes before use.

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