Tetric® PowerFlow

[en] Instructions for Use

Polymer-based dental restorative material
(intraoral light-curing)

Rx ONLY

CE 0123

Date information prepared:
2022-12-20 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Direct restorations of posterior teeth

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

Tetric® PowerFlow is a flowable, light-curing, radiopaque composite (225% Al) for the direct restorative treatment of posterior teeth (according to ISO 4049:2019 Type 1, Class 2, Group 1).

Tetric PowerFlow is also suitable for restoring occlusal surfaces. 100% aluminium has a radiopacity equivalent to that of dentin and 200% aluminium is equivalent to enamel. As its opacity increases during polymerization, Tetric PowerFlow is also suitable for discoloured tooth structure. Due to esthetic reasons, it is applied as initial layer in increments of up to 4 mm in Class I and II restorations.

Tetric PowerFlow cures with light in the wavelength range of 400–500 nm.

Indications

Missing tooth structure in posterior teeth (Classes I and II)

Types of restorations

Conventional application (light intensity $\leq 2,000$ mW/cm²)

- As initial layer / first increment in Class I and II composite restorations in permanent posterior teeth
- Restorations of deciduous teeth

Light-curing using the 3sCure mode of the Bluephase® PowerCure (light intensity 3,000 mW/cm²)

- As initial layer / first increment in Class I and II composite restorations in permanent posterior teeth when light-cured from the occlusal aspect

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working procedures cannot be applied.
- The product is not suitable for Class V restorations (due to its viscosity).

- Due to esthetic reasons, the product is not suitable for restoring occlusal surfaces.
- Avoid direct exposure of the gingiva, mucous membrane or skin to curing light when using the 3sCure mode.
- Tetric PowerFlow can be used in combination with Tetric® PowerFill, Tetric® Prime, Tetric® EvoFlow and Tetric® EvoCeram. Use the 3sCure curing mode of the Bluephase PowerCure only for the Tetric PowerFill or Tetric PowerFlow materials.
- Apply the product at ambient temperature. Cool temperatures render the material difficult to extrude.
- The use of different application tips can render the material difficult to extrude.

Processing restrictions

The product is not intended to be reprocessed or reused. If Tetric PowerFlow is applied directly in the mouth of the patient, the Cavifil and the application tip of the syringe should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).

Side effects

In rare cases, components of Tetric PowerFlow may lead to sensitization. The product must not be used in such cases. In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions

Phenolic substances such as materials containing eugenol/clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with Tetric PowerFlow must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Clinical benefit

Reconstruction of chewing function

Composition

Barium glass, Bis-EMA, copolymer, aromatic methacrylate, Bis-GMA, ytterbium trifluoride, DCP

Total content of inorganic fillers: 46.4 vol%.

Particle size of inorganic fillers: between 0.15 µm and 15.5 µm.

2 Usage

I. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide (e.g. the Tetric PowerFill shade guide). The shade of the shade tab corresponds with the shade of the composite after the change in opacity has taken place, i.e. after polymerization.

II. Isolation

Adequate relative or absolute isolation is required.

III. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp, internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel the enamel margins in the anterior region. Lightly round any sharp enamel edges in the posterior region (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

IV. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. In very deep cavities only, areas close to the pulp should be selectively coated with a calcium hydroxide liner. Subsequently cover with a pressure resistant cement. Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

V. Placement of matrix / interdental wedge

Use a wrap-around matrix or a sectional matrix band for cavities affecting the proximal area and secure it with wedges.

VI. Conditioning / Application of the bonding agent

- Apply the bonding agent to the tooth structure according to the instructions for use of the manufacturer of the product used. Ivoclar recommends using a universal adhesive.

VII. Application of Tetric PowerFlow

- For optimum results, apply Tetric PowerFlow in layers of max. 4 mm thickness using the Cavifil Injector or a comparable applicator.
- Keep the tip of the cannula immersed in the material during application to avoid the formation of bubbles.
-  Tetric PowerFlow must be covered with a layer of methacrylate-based universal or posterior composite. Processing and finishing should be conducted according to the respective instructions for use.
- Proximal contacts must be established by means of a matrix system. The filling material does not exert any pressure on the matrix band. The matrix band can be shaped with a suitable hand instrument (e.g. a balltype condenser) before and/or during light-curing.
- Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
- For the recommendations regarding exposure time (Exposure time) per increment and light intensity (Light intensity) see table 1 (Table 1).
-  **The instructions for use of the curing light must be observed.**
- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again.
- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.

VIII. Finishing / Checking the occlusion / Polishing

Remove any excess material with tungsten carbide or diamond finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Application notes

- Avoid direct exposure of the gingiva, mucous membrane or skin to the light emitted by the curing light.
- The 3sCure mode must not be used in case of caries profunda and very deep cavities.
- The recommended increment thickness is based on hardness profile measurements.
- Syringes or Cavifils must not be disinfected with oxidizing disinfection agents.
- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL008JQ

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Unpolymerized Tetric PowerFlow should not come in contact with skin/mucous membrane or eyes. Unpolymerized Tetric PowerFlow may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Loss of the filling may cause accidental ingestion of material and may make additional dental treatments mandatory.
- Inclusion of air bubbles during the placement of the filling may make additional treatment measures that involve loss of tooth structure mandatory.
- Postoperative sensitivity
- Fracture of the filling or chipping

4 Shelf life and storage

- Storage temperature: 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Expiry date: see information on Cavifils, syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

Table 1

| Light intensity | Exposure time |
|----------------------------------|---------------|
| 500 – 900 mW/cm ² | 20 s |
| 900 – 1,400 mW/cm ² | 10 s |
| 1,800 – 2,200 mW/cm ² | 5 s |
| 2,700 – 3,300 mW/cm ² | 3 s |

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