VOCO IonoStar® Molar

EN Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 9917-1

Product description:

IonoStar Molar is a radiopaque glass ionomer bulk filling material in user-friendly capsules for direct application. In addition to high compressive strength, it boasts good adhesion to enamel and dentine and releases fluoride ions. The usable quantity of cement per capsule is at least 0,15 ml.

Shades: A1, A2, A3

Indications:

- Restorations of non occlusion-bearing class I cavities
- Semi-permanent restorations of class I and II cavities
- Restorations of cervical lesions, class V cavities, root caries
- Restorations of class III cavities
- Restoration of deciduous teeth
- Base/liner
- Core build-up
- Temporary restorations

Contraindications:

IonoStar Molar contains polyacrylic acid. IonoStar Molar should therefore not be used in patients with a known hypersensitivity (allergy) to this ingredient.

Patient target group:

IonoStar Molar is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

Use:

IonoStar Molar should only be used by a professionally trained dental practitioner.

Preparation:

Prepare the cavity using a minimally invasive approach (removing only carious tooth substance). There is no need for undercuts. The layer thickness of the restorations must not be less than 1 mm. Avoid thin tapering margins. Clean the cavity, rinse with water and dab or blow dry. Do not overdry the dentine. It should have a moist, shiny appearance. Place a matrix, if necessary. If the residual dentine thickness is insufficient, the area close to the pulp should be covered with a suitable calcium hydroxide product.

Activating the application capsule and mixing the material:

Only extract the application capsule from the tubular bag immediately before use. During activation, the liquid from the plunger must be pressed into the inside of the capsule as fully as possible. To achieve this, the end of application capsule where the plunger of liquid protrudes is slowly pressed down onto a hard surface (e.g. a table) as far as the plunger will move (see fig. 1). The capsule must only be activated directly before use and loaded immediately into a high-frequency mixer with at least 4,000 oscillations per minute. Mix the activated application capsule for 10 - 15 s (mixing for longer reduces the working and setting times), remove it from the mixer without delay and insert it in the holder of an application device (e.g. AC Applicator, see corresponding instructions for use).

Application:

Carefully open the outlet by lifting the cannula until you feel a stop (see fig. 2), and apply the material directly into the cavity within 15 s of the end of mixing by compressing the lever on the applicator slowly and evenly. The cannula can be aligned as required by turning the whole application capsule inside the holder. Working time and consistency may be influenced by the respective type of capsule mixer and mixing time.


Fig. 1: Activating the capsule Fig. 2: Opening the capsule

Following mixing, IonoStar Molar can be worked for at least 1.5 minutes (at room temperature). Prevent the material from being contaminated with moisture and from drying out. If it is not possible to keep the working area dry, a protective varnish should be applied after the contours/modelling have been finished (see corresponding instructions for use). The restoration can be finished, e.g. with a diamond bur or a polisher, 5 - 6 minutes after the start of mixing at the earliest. A protective varnish should then be applied again (see corresponding instructions for use).

Mixing time (23 °C)	10 - 15 s
End of mixing to application (23 °C)	15 s
Working time (23 °C)	1.5 min
Setting time (37 °C)	3 - 4 min
Total setting time	5 - 6 min

Warnings, precautionary measures:

- Do not apply directly to dentine near the pulp or to exposed pulp.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- Rinse the instruments used immediately with cold water before residues of material cure.
- Ensure that the application capsule is mixed immediately after activation.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Fluorosilicate glass, polyacrylic acid, water, tartaric acid, fumed silica, pigments

Storage instructions and application method:

Store IonoStar Molar in a dry place at temperatures of 15 °C - 23 °C. Only remove it from the foil immediately before use. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of IonoStar Molar must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of IonoStar Molar can be found in the European database on medical devices (EUDAMED - <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 BB2520 E1 0223 99 © by VOCO